


Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method

DETECTION OF VISA, GISA AND hVISA USING THE NEO-SENSITABS PREDIFFUSION METHOD

1. One Neo-Sensitabs of each Vancomycin 30 µg and Teicoplanin 30 µg are placed on an uninoculated plate containing the susceptibility test medium Mueller-Hinton Agar.
2. After 2 hours at room temperature, the tablet (disc) is removed (by knocking the plate against the table), but prior to this a short name (VAN or TEI) is written on the back of the plate in order to make it possible to identify the antimicrobial.
3. Now the plate is maintained at room temperature for a further 18 to 22 hours (overnight).
4. The plate is now inoculated with the strain to be tested using a McFarland 0.5 inoculum. Additional antimicrobial discs (Neo-Sensitabs) may be added now, using a dispenser and thereafter the plate is incubated overnight at 35 degrees.
5. The zones of inhibition are then measured and compared with the corresponding zone breakpoints



S. aureus CB 182, Susceptible strain. V = Vancomycin, T = Teicoplanin, D = Daptomycin

Please note:

In the laboratory the prediffusion plate can be prepared the day before it is inoculated to avoid loss of time and results are available within 24 hours. Surplus of prediffused plates may be kept in the refrigerator for another 24 hours.

Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method

Interpretation

DETECTION of VISA, GISA, hVISA and hGISA

VISA, GISA and hVISA strains will show the following zones of inhibition using Mueller Hinton Agar and McFarland 0.5 inoculum and the prediffusion method:

hVISA, hGISA	VISA, GISA
Teicoplanin 30 ug inhibition zone < 20 mm OR Vancomycin 30 ug inhibition zone ≤ 22 mm	Teicoplanin 30 ug inhibition zone < 20 mm AND Vancomycin 30 ug inhibition zone ≤ 22 mm

Please notice that Teicoplanin in general is the most sensitive drug to detect these isolates. The current MIC methods and automatic systems are unable to detect hVISA strains, because they use too small inocula and consequently cannot detect heteroresistant isolates.

In the laboratory it may be useful to test vancomycin, teicoplanin and daptomycin together in order to be able to detect hVISA/VISA strains first and thereafter in order to find the best drug for treatment of MRSA and hVISA/VISA infections.

A 9 cm Mueller-Hinton Agar plate will be adequate for testing the 3 antimicrobials by the prediffusion method.



S. aureus ATCC 700698 Daptomycin, susceptible hVISA strain.

Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method

DETECTION OF COAGULASE NEGATIVE GLYCOPEPTIDE INTERMEDIATE STAPHYLOCOCCI

The same procedure as above is used for Teicoplanin 30 ug. Isolates showing zones of inhibition < 20 mm with Teicoplanin 30 ug (prediffusion method) should be reported as resistant to teicoplanin and possibly heteroresistant to vancomycin.




S. aureus AT 403, GISA strain, Daptomycin susceptible.

Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar


Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method



Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Prediffusion method (2+18 or 2+22 hours) for antimicrobials diffusing poorly on agar

Detection of VISA, GISA and hVISA using the Neo-Sensitabs prediffusion method



References

- 1) Nielsen SV, Casals JB: Detection of decreased susceptibility to glycopeptides in *S. aureus* using tablet(disc) prediffusion. 15th ECCMID, april 2005.
- 2) Ferreira Nunes AP et al: Heterogeneous resistance to vancomycin in *S. epidermidis*, *S. haemolyticus* and *S. warneri* clinical strains: characterisation of glycopeptide susceptibility profiles and cell wall thickening. *Int J Antimicrob Agents* **27**,307-315,2006.
- 3) Liu C et al: *Staphylococcus aureus* with heterogeneous resistance to vancomycin: epidemiology, clinical significance and critical assessment of diagnostic methods. *Antimicrob Agents Chemother.* **47**,3040-3045,2003
- 4) Sandberg A et al: First report of heterogeneous vancomycin intermediate *S. aureus* (hVISA, VISA) in Denmark. Poster AR-07.ICAAC 2004.
- 5) Jae-Hoon Song et al: Emergence in Asian countries of *S. aureus* with reduced susceptibility to vancomycin. *Antimicrob Agents Chemother* **48**,4926-28,2004.
- 6) Denis O et al: Emergence of vancomycin-intermediate *S. aureus* in a Belgian hospital: microbiological and clinical features. *J Antimicrob Chemother.* **50**,383-391,2002.
- 7) Lodise: Relationship between vancomycin MIC and failure among patients with MRSA bacteremia treated with vancomycin. *Antimicrob Ag. Chemother* **52**,3315-3320,2008.
- 8) In-Gyu Bae et al: hVISA phenotype in bloodstream MRSA isolates from an international cohort of patients with infective endocarditis: prevalence, genotype and clinical significance. *J. Infect Dis* **200**,1355-1366,2009.
- 9) Musta AC et al: Vancomycin MIC plus heteroresistance and outcome of MRSA bacteremia: trends over 11 years. *J. Clin Microbiol* **47**,1640-44,2009.
- 10) Lulitanond A et al: The first VISA strains isolated from patients in Thailand. *J. Clin. Microbiol* **47**,2311-16, 2009.
- 11) Sader HS et al: Occurrence of vancomycin-tolerant and hVISA strains among *S. aureus* causing bloodstream infections in 9 USA hospitals. *J. Antimicrob Chemother* **64**,1024-28,2009.
- 12) Pillai SK et al: Development of reduced vancomycin susceptibility in MSSA. *Clin Infect Dis* **49**,1169-1174,2009.
- 13) Katayama Y. et al: Selection of heterogeneous intermediate *S. aureus* by imipenem. *Antimicrob. Ag. Chemother*, **53**, 3190-96, 2009.
- 14) Fusco D. N. et al: Clinical failure of vancomycin in a dialysis patient with methicillin susceptible, vancomycin-heteroresistant *S. aureus*. *Diagn. Microbiol. Infect. Dis.* **65**, 180-183, 2009.
- 15) Perichon B. et al: Van-A type vancomycin resistant *S. aureus*. *Antimicrob. Ag. Chemother* **53**, 4580-87, 2009.