

Carbapenemases

Carbapenemases are beta-lactamases that significantly hydrolyze at least imipenem and/or meropenem. Carbapenemases involved in acquired resistance are of Ambler classes A, B and D. They may be plasmid or chromosomally encoded.

Because several of these carbapenemases confer only reduced susceptibility to carbapenems in Enterobacteriaceae, they may remain underestimated, because they are not detected in the laboratory. Acquired carbapenemases are increasingly reported worldwide and consequently it is important to be able to detect them in the laboratory.

For many isolates with carbapenemases, the MICs of carbapenems are around the susceptible breakpoint making resistance difficult to detect - particularly with automated systems. Therefore, special zone breakpoints are needed in first line screening

Enterobacteriaceae with reduced susceptibility to Imipenem 10 µg (zone < 23 mm or MIC > 1 µg/ml) or Meropenem 10 µg (inhibition zone < 25 mm or MIC > 0.5 µg/ml) or Ertapenem (zone ≤ 22 mm) on Mueller-Hinton Agar with McFarland 0.5 inoculum, should be suspected of possessing carbapenemases

P. aeruginosa with inhibition zones Imipenem 10 µg (< 22 mm) or Meropenem 10 µg (< 26 mm) should be suspected of possessing carbapenemase. Most isolates with KPC and GES enzymes are highly resistant to Ceftazidime. Ertapenem Neo-Sensitabs is the most sensitive indicator for possible carbapenemase, but in approximately 20% of cases other resistance mechanisms are involved (confirmation of carbapenemase with Modified Hodge Test is necessary). It is important to recognize small resistant colonies growing inside the Ertapenem disk zone.

Carbapenemases classification (1)

| Ambler classification | Enzymes | MICs µg / ml | | | | Inhibited by | | |
|------------------------------|----------------|----------------|-------------|-------------|-----------|--------------|------------|--------------|
| | | 3rd gen cepha | AZT | IMP | MRP | CLAV | EDTA | Boronic acid |
| A | NmcA | S | 4 | ≥ 16 | 2-8 | ± wk | no | yes |
| | Sme-1 to Sme-3 | S | 4-64 | ≥ 16 | 0.25-8 | ± wk | no | yes |
| | IMI-1 to IMI-2 | S | S | ≥ 64 | 4-32 | + | no | yes |
| | KPC-1 to KPC-4 | ≥ 32 | ≥ 64 | 4→16 | 4→16 | + or wk | no | yes |
| | GES-2 to GES-5 | ≥ 32 | 16→R | 0.25→16 | 0.5-16 | + or 0 | no | yes |
| B Metallo-beta-lactamases | IMP 1-16 | ≥ 32 | S→R | 0.5-128 | 0.25→R | no | yes | no |
| | VIM 1-12 | ≥ 64 | S→R | 1→R | 0.5→R | no | yes | no |
| | SPM-1 | ≥ 256 | 4 | R | R | no | yes | no |
| | GIM-1 SIM-1 | 16-32 ≥ 256 | 8-16 128 | > 8 8-16 | > 8 16 | no no | yes yes | no no |
| D Oxacillinases | OXA 23-27 | > 256 | > 256 | 4-64 | 4-128 | ± wk | no | no |
| | OXA 40-48 | S→R | S→R | 2-64 | 0.25-64 | wk | no | no |
| | OXA 54-55 | S | S | 4 | 0.25 | wk | no | no |
| | OXA-60 | S | R | 0.5 | 2 | no | no | no |
| | OXA-58 | 4-128 | ≥ 32 | 3-32 | 2→64 | no | no | no |

wk = weak

References:

- 1) Nordmann P et al: Emerging carbapenemases in gram-negative aerobes. Clin Microbiol. Infect **8**, 321-331, 2002.

Detection of acquired carbapenemases Ambler classes A and D

Class A carbapenemases are penicillinases with greater activity against imipenem than meropenem and they also give resistance to penicillins, cephalosporins and aztreonam.

Boronic acid in an inhibitor of class A carbapenemases and consequently synergy with meropenem or imipenem, is the best method to detect these enzymes (26,27,28,22).

Clavulanate is an inhibitor of class A carbapenemases and therefore synergy with imipenem may be useful to detect these enzymes (1,2,3,4,5).

The KPC family of enzymes confer greater resistance to third gen cephalosporins than to carbapenems (3,5).

KPC possessing *Enterobacter* spp. and *K. pneumoniae* were reported as falsely susceptible to carbapenems using automated systems (Vitek). MIC microdilution using standard inocula of 10⁴ or 10⁵ CFU/ml did not detect carbapenem resistance, while diffusion methods (E-test) using inocula of 10⁸ CFU/ml detected resistance (5,7,12,18).

K. pneumoniae intermediate or resistant to ertapenem or meropenem should be considered resistant to all carbapenems (7). KPC possessing E.coli was identified in nine patients in New York. Three of the isolates possessed also ESBL: CTX M15 (19).

Pasteran et al (20) found that Boronic acid disks could be used to detect carbapenemases of type 2f (Class A) in Enterobacteriaceae. Class A producing strains showed synergy between Imipenem and Boronic acid disks (distance from edge to edge 6 mm). Strains showing zones of inhibition ≤21 mm with Imipenem 10 µg disks were screened with this test.

Carbapenemase IMI-2 is the first inducible and plasmid-encoded carbapenemase.

Please note that KPC detection may require screening multiple colonies, because carbapenemase susceptible strains may co-exist with resistant (21).

Class D carbapenemases correspond to the enzymes classified as OXA-types (oxacillinase activity). They hydrolyze imipenem and meropenem weakly and do not hydrolyze third gen cephalosporins and aztreonam (although MICs against the later drugs are often increased due to the presence of other beta-lactamases).

Clavulanate is a progressive inhibitor of most OXA carbapenemases, but not all. The synergy test (clavulanate and imipenem) may have value for the detection of these enzymes.

Clinicians should be aware of the potential for clinical failure (Class D, OXA-55 carbapenemase) when imipenem is used for treatment of serious infections caused by *S. algae* (9).

Fernandez et al report heteroresistance carbapenems in *Ac. baumannii* associated with Imipenems MIC 4-16 µg and the presence of OXA-58 (17).

Yilmaz et al (16) report oxacillinases (OXA-48) in 21 Enterobacteriaceae, mainly *K.pneumoniae*, but also in *E.coli* and *Enterobacter cloacae/aerogenes* in Turkey, and warns that oxacillinases (carbapenemases) are spreading in Enterobacteriaceae.

Castanheira et al (25) report the clonal dissemination of OXA-24 and OXA-58 producing *A. baumannii* in Houston, Texas.

Detection of resistance mechanisms using Neo-Sensitabs™ and Diatabs™

Detection of beta lactamases

Carbapenemases. KPC+ MBL Confirm ID kit

Acquired carbapenemases Ambler class A and D

| Ambler class. | Enzymes | MICs µg/ml | | | | | IMIPENEM + CLAV (synergy) | Mero-penem+ Boronic synergy | Organisms | Genetic Location |
|---------------|---|------------|---------------|-------|-------|---------|---------------------------|------------------------------------|---|------------------------|
| | | PIPER | 3rd gen cepha | AZTRM | IMI | MEROP | | | | |
| A | NmcA | S | S (0.25-2) | 4 | ≥ 16 | 2-8 | ± wk | yes | <i>E. cloacae</i> | Chromosomal |
| | Sme-1 to Sme-3 | S | 0.25-0.5 | 4-64 | ≥ 26 | 0.25-8 | ± wk | yes | <i>S.marcescens</i> | Chromosomal |
| (11) | IMI-1 | > 256 | 1-2 | 8 | > 64 | 4 | + | yes | <i>E. cloacae</i> | Chromosomal |
| | IMI-2 | 16→128 | 0.1-2 | 4-8 | 64 | 4-32 | + | yes | <i>E. asburiae</i> | Plasmid |
| | KPC-1 | > 128 | ≥ 32 | > 64 | 16 | 16 | + | | <i>K. pneumoniae</i> | Plasmid |
| | KPC-2 | ≥ 64 | ≥ 8 | > 16 | 8→16 | ≥ 16 | + | | <i>K. pneumoniae /oxytoca.</i> <i>Raoultella(39)</i> <i>Salmonella</i> <i>Enterobacter</i> | Plasmid |
| | KPC-2 | . | ≥ 256 | > 256 | 256 | 256 | (+) | yes | <i>P. aeruginosa</i> | Plasmid/Chrom. |
| | KPC-3 | 256 | 256 | > 256 | > 4 | > 4 | (+ +) | | <i>K.pneumoniae</i> <i>Enterobacter</i> | Plasmid |
| | KPC-4 | | | | > 16 | > 16 | (+) | | <i>E. coli, P. mirabilis, Citrobacter</i> | Enterobacter |
| | KPC-5/6 | | | | | | | | | |
| | GES-2 | 128 | ≥ 32 | 16 | 4→16 | 4-16 | + | yes | <i>E. cloacae</i> <i>P. aeruginosa</i> | Plasmid, integron |
| | GES-3 | 128 | 64→256 | 64 | 0.25 | 0.5 | + | yes | <i>K. pneumoniae</i> | Plasmid |
| | GES-4 | 128 | R | R | 8 | 8 | (+) | yes | <i>K. pneumoniae</i> | CEFOX R |
| GES-5 | R | R | R | 8-32 | 8-32 | + | yes | <i>P. aeruginosa, K.pneumoniae</i> | Integron | |
| GES-11 | R | > 256 | > 256 | 4 | 8 | | yes | <i>A. baumannii</i> | | |
| D | OXA-23 to OXA-27 | > 256 | > 256 | > 256 | 4-64 | 4→128 | ± wk | no | <i>A. baumannii</i> | Chromosomal ± integron |
| | OXA-40 | R | 4→128 | 4→128 | > 32 | ≥ 32 | wk | no | <i>Ac.haemolyticus</i> <i>K. pneumoniae</i> <i>P. aeruginosa</i> | Plasmid |
| | OXA-48 | 8→R | S→R | S→R | 2→64 | 0.25→64 | wk | no | <i>K. pneumoniae</i> <i>E.coli</i> | Plasmid |
| | OXA-54 | 32 | S | S | 1 | 0.12 | wk | no | <i>Sh. putrefaciens</i> | Not integron |
| | OXA-55 | S | S | S | 1-4 | 0.25 | no wk | no | <i>Sh. algae (9)</i> | Chromosomal |
| | OXA-58 | 256 | 4-128 | ≥ 32 | 2-32 | 2→64 | no | no | <i>A. baumannii</i> | Plasmid |
| | OXA-60 | S | S | R | 0.5 | 2 | no | no | <i>R. pickettii</i> | Chromosomal |
| | OXA-62 | S→R | S→R | S→R | 2→64 | 64→128 | no | no | <i>Pandorea (10)</i> <i>pnomenusa</i> | Chromosomal |
| D (8) | OXA-23, 27, 49 (subgroup 1) | | | >16 | 8→32 | 8→>32 | | no | <i>Ac. baumannii</i> | Plasmid (only 23) |
| | OXA-24, 25, 26, 40 (subgroup 2) | - | > 256 | | > 128 | > 128 | - | no | <i>Ac. baumannii</i> | Chromosomal |
| | OXA-51 + OXA-64-66, 68-71, 78-82-107 OXA-51-like (subgroup 3) | | | | ≥ 1 | ≥ 1 | | no | <i>Ac. baumannii</i> | Chromosomal plasmid |
| | OXA-58 (subgroup 4) | R | R | > 16 | 4/16 | | | no | <i>Ac. baumannii</i> | Plasmid (only 58) |
| | OXA-143 | R | FEP4 | . | 32 | 32 | - | no | <i>Ac. baumannii</i> | Plasmid |

Bold = involved in outbreaks

Procedure for KPC carbapenemases detection (Class A enzymes)

Isolates giving negative metallo-beta-lactamase tests, may produce other carbapenemases. The most current are KPC enzymes isolated from *Enterobacteriaceae* (*K. pneumoniae*, *E. coli*, *Enterobacter spp.*, *P. mirabilis*) particularly *K. pneumoniae*, but also Sme, IMI, GES and Nunc A are found. To detect these strains in rectal swab screening samples, direct plating on McConkey agar in the presence of Ertapenem Neo-Sensitabs and Imipenem Neo-Sensitabs may be useful.

Place one Boronic Acid Diatabs between one Ertapenem and one Imipenem Neo-Sensitabs (distance 6 mm from edge to edge).

Place one Cloxacillin Diatabs between Ertapenem and Imipenem Neo-Sensitabs - (6 mm from edge to edge). Perform Modified Hodge Test (MHT) with Ertapenem and Meropenem Neo-Sensitabs.


Interpretation (Double disk synergy test)

The following results will presumably indicate the presence of a KPC beta-lactamase:

- a) Negative metallo-beta-lactamases tests.
- b) Positive synergy test between Boronic Acid and the carbapenems (one or both).
- c) Negative synergy test between Cloxacillin and the Carbapenems (11)
- d) Positive synergy test between clavulanate (AMC) and carbapenems (one or both). Not always easy to see. Although isolates with ESBL + impermeability may give false positive results.
- e) Positive Modified Hodge Test.
- f) Sme, IMI, GES and Nunc A will show the same results as KPC, but the mentioned enzymes result in smaller zones around Imipenem compared to Ertapenem. With KPC enzymes zones around Imipenem and Ertapenem are similar.

Combined disk test. KPC + MBL Confirm ID kit

Apply Meropenem, Meropenem+DPA, Meropenem+Boronic, Meropenem+Cloxacillin on an inoculated plate.



K.pneumoniae PHA3 CL5761 KPC positive

Interpretation (combined test)

A Meropenem + Boronic inhibition zone ≥ 5 mm then Meropenem, Meropenem+DPA and Meropenem+Cloxacillin indicates a presence of a KPC enzyme (or other class A). Meropenem+Boronic and Meropenem+Cloxacillin inhibition zones ≥ 5 mm, than Meropenem and Meropenem+DPA indicates AmpC hyperproduction + porin loss, or efflux (30).

A Meropenem +DPA inhibition zone ≥ 5 mm than Meropenem, indicates the presence of a metallo- β -lactamases (MBL).

Please note

Test only ertapenem-resistant strains. Ertapenem susceptible strains may provide a false positive result with Boronic Acid.

Conclusion

Reduced susceptibility to ertapenem, synergy between Boronic Acid and the carbapenems, and no synergy between Cloxacillin and the carbapenems is clearly indicative of KPC enzyme being present (or other class A enzymes). Isolates producing high level AmpC + impermeability can be detected by synergy between Cloxacillin and the carbapenems (11). Isolates producing ESBL + impermeability will show synergy between AMC and the carbapenems or cephalosporins.

Procedure for Oxacillinase detection (Class D enzymes)

Strains producing oxacillinases will currently show zones of inhibition < 22 mm with Ertapenem and/or <25 mm with Meropenem Neo-Sensitabs. Most are resistant to Aztreonam. These enzymes are mainly found in *Acinetobacter baumannii* but also in Enterobacteriaceae (*K. pneumoniae*, *Enterobacter*) and *P. aeruginosa* although these are rare.

Interpretation

The following results will presumably indicate the presence of oxacillinases:

- a) Negative metallo-beta-lactamase tests.
- b) Negative synergy test between Boronic acid/Cloxacillin and the carbapenems. (one or both).
- c) Negative (or weak positive) synergy test between clavulanate (AMC) and carbapenems (one or both)
- d) Positive Modified Hodge Test.

Summary detection of carbapenemases

| Meropenem | MRP+DPA | MRP+BOR | MRP+Cloxa |
|----------------------|----------------|----------------|------------------|
| Metallo-β-lactamases | Synergy | No synergy | No synergy |
| KPC | No synergy | Synergy | No synergy |
| AmpC impermeability | No synergy | Synergy | Synergy |
| Oxacillinases | No synergy | No synergy | No synergy |


Modified Hodge Test

Is used to determine if resistance to carbapenems is caused by a carbapenemase. A MH agar plate (or a McConkey plate) is inoculated with the susceptible strain *E. coli* ATCC 25922 (Mc Farland 0.5, diluted 1/10) as for disk diffusion.

When testing Enterobacteriaceae, one Ertapenem Neo-Sensitabs and one Meropenem Neo-Sensitabs are applied onto the plate approx.30 mm apart from each other. For non-fermenters one Imipenem Neo-Sensitabs and one Meropenem Neo-Sensitabs are applied.

A suspension of the microorganism to be tested for carbapenemase is adjusted to Mc Farland 0.5 standard and a loop is used to make a heavy streak passing through the two carbapenem disks. Two more streaks are placed perpendicularly making a cross.

Thereafter incubation for 18-24 hours at 35-37 C. Alteration in the shape (indentation) of the zones of inhibition around the test organism is considered indicative of the presence of a carbapenemase (figure).



K.pneumoniae KPC positive

Limitations:

- Not reliable for detection of SME from *S. marcescens*
- *P. mirabilis* swarming may give lecture problems

References:

- 1) Pottumarthy S et al: NmCA carbapenem-hydrolyzing enzyme in *E. cloacae* in North America. *Emerg Infect Dis* **9**, 999-1002, 2003.
- 2) Yigit H et al: Novel carbapenem-hydrolyzing beta-lactamase KPC-1 from a carbapenem-resistant strain of *K. pneumoniae*. *Antimicrob. Ag. Chemother.* **45**, 1151-61, 2001.
- 3) Smith-Moland E et al: Plasmid-mediated carbapenem-hydrolyzing beta-lactamase KPC-2 in *K. pneumoniae* isolates. *J. Antimicrob. Chemother.* **51**, 711-14, 2003.
- 4) Aubron C et al: Carbapenemase-producing Enterobacteriaceae, U.S. rivers. *Emerg. Infect. Dis.* **11**, 260-4, 2005.

- 5) Bratu S. et al: Detection of KPC carbapenem-hydrolyzing enzymes in Enterobacter spp from Brooklyn, N.Y. Antimicrob. Ag. Chemother. **49**, 776-8, 2005.
- 6) Lopez-Otsoa F et al: Endemic carbapenem resistance associated with OXA-40 carbapenemase among *A. baumannii* isolates from a hospital in Northern Spain. J.Clin.Microbiol **40**, 4741-3, 2002.
- 7) Bratu S. et al: Emergence of KPC-possessing *K. pneumoniae* in Brooklyn, N.Y.: epidemiology and recommendations for detection. Antimicrob. Ag. Chemother., **49**, 3018-20, 2005.
- 8) Brown S., Amyes S.: OXA beta-lactamases in Acinetobacter: the story so far. JAC, **57**, 1-3, 2006.
- 9) Dong-Min K. et al: Treatment failure due to emergence of resistance to carbapenem during therapy for *Shewanella alga bacteremia*. J.C.M., **44**, 1172-4, 2006.
- 10) Schneider I. et al: Novel carbapenem-hydrolyzing oxacillinase OXA-62 from *Pandoraea pnomenusa*. Antimicrob. Ag. Chemother., **50**, 1330-35, 2006.
- 11) Villegas M.V. et al: First identification of *P. aeruginosa* isolates producing a KPC-type carbapenem hydrolyzing β -lactamase. Antimicrob. Ag. Chemother., **51**, 1553-55, 2007.
- 12) Tenover F.C. et al: Carbapenem resistance in *K. pneumoniae* not detected by automated susceptibility testing. Emerg. Inf. Dis., **12**, 1209-1213, 2006.
- 13) Sundin D.R.: Rapid screening for KPCs. Presentation D 1555, ICAAC september 2007, Chicago, USA.
- 14) Queenan A.M. et al: Carbapenemases: the versatile β -lactamases. Clin. Microbiol. Reviews, **20**, 440-458, 2007.
- 15) Lee K. et al: Effect of Oxgall on the Imipenem disk Hodge test in screening MBL- producing gram-negative bacilli. ECCMID, poster P-891, Barcelona 2008.
- 16) Yilmaz M. et al: Warning! Oxallinase-mediated carbapenem-resistance spreading in Enterobacteriaceae. Poster P1517, ECCMID. Barcelona, April 2008.
- 17) Fernandez F. et al: Association between heteroresistance to carbapenems in *A. baumannii*, production of beta-lactamases (AmpC and OXA). XIII Congress SEIMC presentation 12, Madrid, May 2008.
- 18) Marchaim D. et al: Isolation of Imipenem-R Enterobacter spp: emergence of KPC-2 carbapenemase, molecular characterisation, epidemiology and outcomes. Antimicrob. Ag. Chemother **52**, 1413-18, 2008.
- 19) Urban C. et al: Carbapenem-resistant *E. coli* harboring *K. pneumoniae* carbapenemase beta-lactamases associated with long term care facilities. CID **46** e127-e130, 2008.
- 20) Pasteran F. et al: Detection of carbapenemases of group 2f. Promising additive uses of the Boronic acid disk test. Poster 371-20836 XI, Argentinian Congress of Microbiology 2007.
- 21) Kohlerschmidt D.J. et al: *Klebsiella pneumoniae* KPC detection may require screening multiple colonies for accurate carbapenem resistance assessment. ASM, presentation C-028, June 2008.
- 22) Pasteran F. G. et al: *Klebsiella pneumoniae* KPC-2, Buenos Aires, Argentina. Emerg. Infect. Dis. **14**, 1178-1180, 2008.
- 23) Ya-Gang Chen et al: In vivo development of carbapenem resistance in clinical isolates of *E. aerogenes*, producing multiple beta-lactamases. Int. J. Antimicrob. Ag, **32**, 302-307, 2008.
- 24) Cuzon G. et al: Plasmid encoded carbapenem hydrolyzing β -lactamase OXA-48 in an Imipenem susceptible strain from Belgium. Antimicrob. Ag. Chemother, **52**, 3463-4, 2008.
- 25) Castanheira M. et al: Emergence and clonal dissemination of OXA-24 and OXA-48 producing *A. baumannii* strains in Houston, Texas: report from the SENTRY-ASP. J. Clin. Microbiol. **46**, 3179-80, 2008.
- 26) Doi Y. et al: Simple disk based method for detection of *K. pneumoniae* carbapenemase by use of Boronic acid compound. J. Clin Microbiol. **46**, 4076-83. 2008.
- 27) Giske C. et al: Personal communication 2009
- 28) Tsakris A. et al: First occurrence of KPC-2 possessing *K. pneumoniae* in a Greek hospital and recommendation for detection with Boronic acid disk tests. J. Antimicrob. Chemother **62**, 1257-60, 2008.
- 29) Pournaras S. et al: Clonal spread of KPC-2 carbapenemase producing *K. pneumoniae* in Greece. JAC, **64**. 348-52, 2009.
- 30) Smith Moland E. et al: Concerns about KPC screening and confirmatory tests. Presentation D-729, 49th ICAAC 2009.
- 31) Endimiani A. et al: Evaluation of the 2009 CLSI criteria for determining carbapenem susceptibility of *K. pneumoniae*. Presentation D-722, 49th ICAAC, 2009.
- 32) Anderson K. et al: Characterization of Enterobacteriaceae with a false positive Modified Hodge Test. Presentation D-719, 49th ICAAC, 2009.

- 33) Landman D. et al: Contribution of OmpK 36 to carbapenem susceptibility in KPC-producing *K. pneumoniae*. J. Med. Microbiol. **58**, 1303-1308, 2009.
- 34) Kristo et al: Investigation of meropenem heteroresistance in KPC-2 producing *K.pneumoniae*. Presentation E-183. 49th ICAAC. 2009.
- 35) Tsakris A. et al: Evaluation of Boronic acid disk tests for differentiating KPC-possessing *K. pneumoniae* isolates in the clinical laboratory. J. Clin. Microbiol. **47**, 362-67, 2009.
- 36) Pasteran F. et al: Sensitive screening tests for suspected class A carbapenemase production in species of Enterobacteriaceae. J. Clin. Microbiol. **47**, 1631-39, 2009.
- 37) Weisenberg S.A. et al: Clinical outcomes of patients with *K. pneumoniae* KPC producing, after treatment with imipenem or meropenem. Diagn. Micr. Infect. Dis. **64**, 233-35, 2009.
- 38) Nicoletti A.G. et al: *K. pneumoniae* non-susceptible to carbapenems in Brazil. What are the mechanisms behind? Presentation C2-663, 49th ICAAC, 2009.
- 39) Castanheira M et al: First descriptions of KPC in *Raoultella* spp (*R. planticola* and *R. ornitholytica*): Report from the Sentry Antimicrobial Surveillance program. J. Clin. Microbiol. **47**, 4129-30, 2009.

Detection of acquired Metallo-beta-lactamases (MBL)

The worldwide spread of acquired metallo-beta-lactamases (MBL) in gram-negative aerobes is of great concern. MBL production in clinical isolates of key gram-negatives: *P. aeruginosa*, *E. cloacae*, *S. marcescens* and *K. pneumoniae* should be carefully monitored (5).

MLBs are classified into 5 major types: IMP, VIM, SPM, and GIM and SIM type enzymes. In Enterobacteriaceae only IMP and VIM enzymes have been found as yet.

MBLs hydrolyze most beta-lactams (carbapenems and large spectrum cephalosporins), except aztreonam. This phenotype of multiple beta-lactam resistance and aztreonam susceptibility may be helpful for identification of these strains in the laboratory. If the strain is resistant to aztreonam it may be due to additional resistance mechanisms (efflux, other beta-lactamases, ESBL etc.). Their expression is not inducible.

The MBL enzymes are resistant to beta-lactamase inhibitors and susceptible to chelating agents like EDTA (2-MPA) and Dipicolinic acid (DPA).

Early detection of MBL-producing microorganisms is essential to prevent dissemination of these organisms. The enclosed tables, including strains of Enterobacteriaceae and Non-fermenters producing MBLs, show that MBL-producers (particularly in Enterobacteriaceae) may show low MIC values against carbapenems making it difficult for the laboratory to detect MBL-positive isolates.

Suspicious isolates (resistant to ceftazidime showing no synergy between clavulanate and third gen. cephalosporins and possibly showing reduced susceptibility to carbapenems) should be tested for carbapenemase activity using Imipenem, Meropenem and EDTA and Dipicolinic acid tests.

The first metallo-beta-lactamase producing strain of *E. coli* (in Spain) has been detected in Barcelona, using Imipenem+EDTA Neo-Sensitabs and E-test (3,8). The first metallo-beta-lactamase producing strain of *K. pneumoniae* was found in France (4).

MBL- producing gram-negatives have now emerged in Australia (15).The resistance gene bla-IMP4 appears highly mobile, this outbreak involved 5 different gram-negative genera. Diagnostic laboratories in Australia and other countries must be now in high alert, because early detection may limit the wide dispersal of MBL-genes.

Kyegong (27) and Miriagou (28) showed the efficiency of Dipicolinic acid (DPA) to detect metallo-β-lactamases in Enterobacteriaceae and non-fermenters. Miriagou found that the DPA/Imipenem synergy test was positive for all VIM-producing isolates of Klebsiella/Enterobacter and *P. mirabilis*, while EDTA based tests could not identify VIM-producing *P. mirabilis*

Acquired Metallo-beta-lactamases NON-FERMENTERS

| MBL | 3rd gen. Cepha MIC | AZT MIC µg/ml | IMP MIC µg/ml | MRP MIC µg/ml | Microorganisms | Genetic location |
|-----------|--------------------|---------------|---------------|---------------|---|--|
| IMP 1-11 | ≥ 128 | ≤ 8/16 | ≥ 8 | ≥ 8 | <i>Pseudomonas</i> spp. <i>Alcaligenes</i> spp. <i>Acinetobacter baumannii</i> | } Chromosomal plasmid integron integron |
| IMP 12 | ≥ 128 | 32 | 32 | 128 | <i>Pseudomonas putida</i> | |
| IMP 13-16 | ≥ 256 | 4-128 | ≥ 64 | ≥ 64 | <i>Pseudomonas aeruginosa</i> | |
| VIM 1-3 | R | S → R | 2-128 | 1-128 | <i>Achromobacter xylosoxidans</i> <i>Pseudomonas aeruginosa</i> <i>Pseudomonas putida</i> (VIM 2 and 4) | } Chromosomal plasmid integron |
| VIM 4-11 | > 256 | S → R | 32-256 | 32-256 | <i>Acinetobacter baumannii</i> <i>Pseudomonas aeruginosa</i> , <i>A.baumannii</i> | |
| VIM 15-16 | ≥ 64 | 16-32 | >128 | ≥128 | <i>Pseudomonas aeruginosa</i> | integron |
| VIM-18 | R | S → R | R | R | <i>Pseudomonas aeruginosa</i> | integron |
| SPM-1 | ≥ 256 | 4 | R | R | <i>Pseudomonas aeruginosa</i> | Plasmid (not integron) |
| GIM-1 | 16 → 32 | 8-16 | > 8 | > 8 | <i>Pseudomonas aeruginosa</i> | Integron |
| SIM-1 | ≥ 256 | 128 | 8-16 | 16 | <i>Acinetobacter baumannii</i> | Integron |
| IND1-6 | 1-32-128 | 32-128 | 4-32-128 | 4-16-128 | <i>Chryseobact indologenes</i> | Chromosomal (23) |
| AIM-1 | | | | | <i>Pseudomonas aeruginosa</i> | |

MBL are not inhibited by clavulanate, but are inhibited by EDTA or DPA

Acquired Metallo-beta-lactamases ENTEROBACTERIACEAE

| MBL | 3rd gen. Cepha MIC | AZT MIC µg/ml | IMP MIC µg/ml | MRP MIC µg/ml | Microorganisms | Genetic location |
|--------------|--------------------|---------------|---------------|---------------|--|-----------------------------|
| IMP-1 | ≥ 32 | < 0.5 | 2 | 0.5 | <i>E. coli</i> | } Integron } Plasmid |
| IMP-1 | ≥ 32 | 0.5 → R | 4-128 | 4-128 | <i>S. marcescens</i> , <i>K. pneumoniae</i> , <i>K. oxytoca</i> , <i>E. cloacae</i> / <i>E. aerogenes</i> , <i>Cit. freundii</i> , <i>P. rettgeri</i> , <i>M. morgani</i> , <i>Shigella flexneri</i> | |
| IMP-3 | 64 | 0.5 | 1 | . | <i>Citrobacter youngae</i> | |
| IMP-4 | 256 | . | 3 | 6 | <i>E. coli</i> | |
| IMP-6 | > 128 | 0.25 | 2-8 | 64 | <i>Serratia marcescens</i> | |
| IMP-6 | > 128 | 128 | 32 | > 128 | <i>Enterobacter cloacae</i> , <i>Klebsiella pneumoniae</i> , <i>S.marcescens</i> | |
| IMP-8 | R | S → R | 0.5-8 | 0.25-4 | | |
| VIM-1 | R | 8-128 | 8-32 | 2-32 | <i>E. coli</i> , <i>P. mirabilis</i> (integron) <i>C. koseri</i> , <i>K. oxytoca</i> <i>Klebsiella pneumoniae</i> , <i>E. cloacae</i> | |
| VIM-1 | 16-128 | S → R | 1-64 | 1-32 | <i>Citrobacter freundii</i> / <i>E.cloacae</i> | Plasmid |
| VIM-2 | ≥ 32 | S → R | ≥ 1 | 0.5 → > 2 | <i>Serratia marcescens</i> , <i>P. rettgerii</i> | Integron |
| VIM-2 | ≥ 128 | 32 | 16-64 | 8-64 | <i>Klebsiella oxytoca</i> | Plasmid (integron) |
| VIM-2 | 8 | 16 | 4 | 0.1 (S) | <i>K. pneumoniae</i> / <i>E. cloacae</i> | Plasmid |
| VIM-4 | ≥ 32 | 4 → R | 2-4 | 0.5-1 | <i>K. pneumoniae</i> | Plasmid (16) |
| VIM-12 | ≥ 128 | 16 | 8 | 4 | <i>E. coli</i> | Plasmid (22) |
| VIM-12 | > 32 | 1 | 1 | 0.25 | <i>E. coli</i> | Integron |
| VIM-2 + GES7 | | | | | | |
| KHM-1 | R | 0.25 | 2 | 4 | <i>C. freundii</i> | Plasmid |
| NDM-1 | R | | | | <i>K. pneumoniae</i> , <i>E.coli</i> , <i>C.freundii</i> (31,33) | Plasmid |

MBL are not inhibited by clavulanate, but are inhibited by EDTA or DPA.

Procedure for metallo-beta-lactamase (MBL) detection

Some resistance profiles may suggest MBL production, for example:

a) *Pseudomonas aeruginosa*, *Pseudomonas* spp. and *Acinetobacter* spp.

All isolates non-susceptible to carbapenems and resistant to either ticarcillin, ticarcillin+clavulanate or ceftazidime should be tested for MBL production.

b) *Enterobacteriaceae*

For *E. coli*, *Klebsiella* spp., *P. mirabilis*, *Salmonella* spp. and *Shigella* spp.: All carbapenem S-I-R isolates that are resistant to ceftazidime and amoxicillin+clavulanate and are non-susceptible to ceftazidime (inhibition zone < 18 mm) should be tested for MBL production. In all other cases all isolates are non-susceptible to carbapenems (18).

Enterobacteriaceae

Apply one Dipicolinic Acid Diatabs (DPA) on an inoculated Mueller Hinton (MH) plate. Apply one Meropenem Neo-Sensitabs and one Ertapenem Neo-Sensitabs onto the plate on either side of the DPA, 5mm from DPA (edge to edge). Apply Imipenem 10 µg + EDTA (IM10E) on an inoculated MH plate. Apply Imipenem 10 µg Neo-Sensitabs.


Non-fermenters

Apply one DPA Diatabs on the MH plate. Apply one Imipenem Neo-Sensitabs and one Meropenem Neo-Sensitabs at either side of the DPA, 5 mm from the DPA (edge to edge). Apply Imipenem + EDTA (IM10E) on the inoculated MH plate. Apply one Imipenem Neo-Sensitabs.

Interpretation (Double disk synergy test)

The use of two chelating agents EDTA and DPA will enhance the detection of metallo-β-lactamases (MBL) in the clinical laboratory. A key hole or ghost zone between carbapenems (one or more) and Dipicolinic Acid indicates the presence of an MBL.


- An Imipenem + EDTA (10+750 µg) zone 7mm larger than Imipenem 10 µg indicates the presence of a Metallo-Beta-Lactamase.



A key hole or ghost zone between carbapenems (one or more) and Dipicolinic Acid indicates the presence of a MBL.

Combined disk test. KPC and MBL Confirm ID kit

Apply Meropenem, Meropenem+DPA on an inoculated MH plate. Interpretation: A Meropenem+DPA inhibition zone ≥ 5 mm than Meropenem alone indicates the presence of a metallo-beta-lactamase.



P.aeruginosa FN 8173 metallo- β -lactamase positive

Please note:

Most MH agar brands contain physiological levels of Zn⁺⁺ ions and should be used for carbapenem testing. Iso-Sensitest agar has low levels of zinc ions and may give false susceptibility results for carbapenems in the presence of MBL. Strains of *Acinetobacter baumannii* producing certain oxacillinases may give a false positive metallo- β -lactamase test result.

References:

- 1) Walsh T.R. et al: Evaluation of new E-test for detecting metallo-beta-lactamases in routine clinical testing. *J. Clin. Microbiol.*, **40**, 2755-59, 2002.
- 2) Jong D. et al: Imipenem-EDTA disk method for differentiation of metallo-beta-lactamase producing clinical isolates of *Pseudomonas* spp. and *Acinetobacter* spp. *J. Clin. Microbiol.*, **40**, 3798-3801, 2002.
- 3) Larrosa M. N. et al: *E. coli* multi-resistente productora de una metalo-beta-lactamasa (Spanish). SEIMC Congress, Bilbao, 16-19th May, 2004.
- 4) Lartigue MF et al: First detection of a carbapenem-hydrolyzing metalloenzyme in an Enterobacteriaceae isolate in France. *Antimicr. Ag. Chemother.*, **48**, 4929-30, 2004.
- 5) Castanheira M et al: Molecular characterization of a beta-lactamase gene blaGIM-1, a new subclass of metallo-beta-lactamases. *Antimicr. Ag. Chemother.*, **48**, 4654-61, 2004.
- 6) Jing-Jou Yan et al: Comparison of the double-disk, combined disk and E-test methods for detecting metallo-beta-lactamases in gram-negative bacilli. *Diagn. Microbiol. Infect. Dis.*, **49**, 5-11, 2004.

- 7) Walsh T.R. et al: Metallo-beta-lactamases: The quite before the storm? Clin. Microbiol. Reviews, **18**, 306-325, 2005.
- 8) Tortola M.T. et al: First detection of a carbapenem-hydrolyzing metalloenzyme in two Enterobacteriaceae isolates in Spain. Antimicrob. Ag. Chemother., **49**, 3492-4, 2005.
- 9) Chu Y.W. et al: EDTA susceptibility leading to false detection of metallo-beta-lactamase in *P. aeruginosa* by E-test and an Imipenem-EDTA disk method. Intl. J. Antimicrob. Agents, **26**, 340-341, 2005.
- 10) Lincopan N. et al: First isolation of metallo-β-lactamase producing multiresistant *K. pneumoniae* from a patient in Brazil. J. Clin. Microbiol., **43**, 516-9, 2005.
- 11) Pitout J.D.D. et al: Detection of *P. aeruginosa* producing metallo-β-lactamases in a large centralized laboratory. J. Clin. Microbiol., **43**, 3129-35, 2005.
- 12) Marqué S. et al: Regional occurrence of plasmid-mediated carbapenem-hydrolyzing oxacillinase OXA-58 in *Acinetobacter* spp. in Europe. J.Clin. Microbiol., **43**, 4885-8, 2005.
- 13) Choi Y.S. et al: Evaluation of Imipenem Disk Hodge Test and double disk synergy tests to detect SIM-1 type metallo-β-lactamase-producing *Acinetobacter* isolates. Presentation D-1705, 45th ICAAC, 2005.
- 14) Lee Kyungwon et al: Novel acquired metallo-β-lactamase gene blaSIM-1 in class 1 integron from *Ac. baumannii* clinical isolates from Korea. Antimicrob. Ag. Chemother., **49**, 4485-91, 2005.
- 15) Peleg A.Y.: Dissemination of the metallo-β-lactamase gene blaIMP-4 among gram-negative pathogens in a clinical setting in Australia. CID, **41**, 1549-56, 2005.
- 16) Pournaras et al: VIM-12 a novel plasmid-mediated metallo-β-lactamase from *K. pneumoniae* that resembles a VIM-1/VIM-2 hybrid. Antimicrob. Ag. Chemother., **49**, 5153-6, 2005.
- 17) Berges L. et al: Prospective evaluation of Imipenem/EDTA combined disk and E-test for detection of MBL producing *P. aeruginosa*. Clin. Microbiol. Infect., **12**, Suppl. 4, P1451, 2006.
- 18) Gornaglia G. et al: Metallo-β-lactamases as emerging resistance determinants in Gram-negative pathogens; open issues. Int. J. Antimicrob. Ag., **29**, 380-88, 2007.
- 19) Franklin C. et al: Phenotypic detection of carbapenem-susceptible metallo-β-lactamase-producing Gram-negative bacilli in the clinical laboratory. J. Clin. Microbiol., **44**, 3139-3144, 2006.
- 20) Moland E.S. et al: Prevalence of newer β-lactamases in gram-negative clinical isolates collected in the U.S. from 2001 to 2002. J. Clin. Microbiol., **44**, 3318-3324, 2006.
- 21) Deshpande L.M. et al: Emergence of serine carbapenemases (KPC and SME) among clinical strains of Enterobacteriaceae in the U.S. Medical Centers: Report from the MYSTIC Program (1999-2005).
- 22) Iconomidis A. et al: First occurrence of an *E. coli* clinical isolate producing the VIM1/VIM2 hybrid metallo-β-lactamase VIM12. Antimicrob. Ag. Chemother., **51**, 3038-39, 2007.
- 23) Perilli M. et al: Identification and characterisation of a new metallo-β-lactamase from a clinical isolate of *Chryseobacterium indologenes*. Antimicrob. Ag. Chemother., **51**, 2988-2990, 2007.
- 24) Anderson K.F. et al: Evaluation of methods to identify the *K. pneumoniae* carbapenemase in Enterobacteriaceae. J. Clin. Microbiol., **45**, 2723-25, 2007.
- 25) Kim Soo-Young et al: Convenient test using a combination of chelating agents for the detection of metallo-beta-lactamases in the clinical laboratory. J. Clin. Microbiol., **45**, 2798-2801, 2007.
- 26) Giakkoupi P. et al: Supplementation of growth media with Zn²⁺ facilitates detection of VIM-2 producing *P. aeruginosa*. JCM **46**, 1568-9, 2008.
- 27) Kyeong S.S. et al: Dipicolinic acid and based disk methods for detection of metallo-β-lactamase producing *Pseudomonas* spp and *Acinetobacter* spp. Diagn. Microbiol. Infect Dis **62**, 102-105, 2008.
- 28) Miriagou V et al: Improved detection of VIM-1 production in *Proteus mirabilis* by an Imipenem-dipicolinic acid double disk synergy test. In press.
- 29) Multiresistant hospital bacteria linked to India and Pakistan. HPS Weekly Report News, July 8th, 2009.
- 30) Carvalho-Assef APD et al: First report of the carbapenem-resistant IMP-16 producing *Pseudomonas stutzeri*. Presentation C2-634, 49th ICAAC, 2009.
- 31) Pasteran F. et al: Spread of multiple clones of metallo-β-lactamases (MBL) producing *Pseudomonas putida*. Presentation C2-636, 49th ICAAC, 2009.
- 32) Treviño M. et al: Carbapenem-resistant *E. cloacae* and the emergence of metallo-β-lactamases producing strains in a third level hospital (Santiago de Compostela, Spain). Eur. J. Clin. Microbiol. Infect. Dis **28**, 1253-58, 2009.

- 33) Dongeun Yong et al: Charatcterization of a new metallo- β -lactamase gene NDM-1 and a novel erythromycin esterase gene carried on a unique genetic structure in *K. pneumoniae* sequence type 14 from India. *Antimicrob. Ag. Chemother* **53**, 5046-5054, 2009.